

Sample Collection Equipment List (BLOOD ONLY)


Below are the details of equipment required to conduct blood sample collections. Please refer to these lists to ensure you take all the equipment required to conduct each test.

1. Blood Collection (Capillary Finger Prick) for PEth test

1. SOPDRG009S Collection Instructions
2. Cansford Blood Kit containing the following: (ALWAYS take spares)
 - Sample envelope
 - Hemaxis blood collection card **Check expiry date**
 - Alcohol swab
 - Gauze
 - Lancet
 - Plaster
 - Special delivery return envelope
 - Tamper seal
 - Gloves
3. Hair/ Blood Collection and Consent form
4. Spare Lancets
5. Sharps Bin (this should be fully filled before returning to CHG)

2. Venous Blood Collection for Synlab testing

1. SOPDRG012S Collection Instructions
2. Synlab Blood Kit containing the following: (ALWAYS take spares)
 - Synlab return envelope
 - Synlab test order form (test has been selected using attached label)
 - 21G green needle
 - Protective clam shell for blood tubes
 - EDTA vacutainer tube (purple) and SST vacutainer tube (yellow) **Check expiry date**
 - Plastic vacutainer holder
 - Alcohol Wipe
 - Cotton wool
 - Plaster
3. Hair/ Blood Collection and Consent form
4. Gloves
5. Tourniquet
6. Kidney bowl
7. Sharps Bin (this should be fully filled before returning to CHG)
8. CHG prepaid 1st Class return envelope
9. 22G black needle & 23G butterfly needle for small / difficult veins

 CRYSTAL Health Group	Document name	Version
	Blood Sample Collection Equipment List	v03.04.25